



Tote Agitation Instructions for

SOIL AMENDMENT | ORGANIC SOIL AMENDMENT | FX FOLIAR TREATMENT

IMPORTANT TO READ

If a persistent layer remains on the inside of the tote after agitation, there is no cause for concern. PhycoTerra® is an inert microalgae product, and this residue will not affect performance in the tank or in the field.

Easy Steps

AT INITIAL USE

For agitation of 250-gallon tote, ensure that mixing shaft does not exceed 32" in and mixer power is at a minimum of 0.5 hp. Clean the shaft and blade prior to use. For minimal compaction, agitate for 20 minutes. For observed compaction agitate a minimum of 30 minutes.

NOTE

If entire tote is not used upon first use, continue with steps 1-3, below, before each daily use until the tote is exhausted.

STEP 1

Before each use, agitate contents. Some settling within the tote is normal.

STEP 2

Agitate using a mixing shaft or recirculate with a circulation pump. This will help reconstitute the product into solutions.

STEP 3

Empty tote contents into tank, and repeat steps 1 - 3 until product is exhausted.